

NeuroSense Tremor Classifier — Full Plan

Stage 1 — Data Loading & Subject Selection

- Load `imu_sdata.pickle` and `subject_metadata.pickle`
 - Split subjects by `healthstatus_id`: PD patients (0) vs Healthy Controls (1 and 2)
 - For each subject, flatten all sessions into a list of 5-second windows with metadata attached:
 - subject ID
 - health status
 - session index
 - window index within session
 - UPDRS tremor subscore if available (from `med_eval_1` / `med_eval_2`)
-

Stage 2 — Feature Extraction Per Window

Each window is shape `(500, 3)` — 500 samples × 3 axes at 100Hz.

Step 1 — Preprocessing

- Compute signal magnitude: `mag = sqrt(x2 + y2 + z2)` — collapses 3 axes into one signal, more robust to phone orientation
- Apply a high-pass filter (>1 Hz) to remove gravity component

Step 2 — Spectral Analysis

- Run FFT on the magnitude signal
- Compute Power Spectral Density (PSD)

Step 3 — Extract 3 Features

- **Tremor power** — sum of PSD in the 4–6 Hz band
 - **Dominant frequency** — frequency with peak power within 3–7 Hz (slightly wider to catch edge cases)
 - **Tremor regularity** — sharpness of the spectral peak, computed as ratio of peak power to mean power in the 3–7 Hz band. High ratio = narrow sharp peak = sustained rhythmic tremor
-

Stage 3 — Threshold Calibration from SData

This is what makes the classifier data-driven rather than arbitrary.

- Compute tremor power for every window across **all PD subjects only**
 - Build the distribution
 - Set thresholds:
 - **LOW_THRESH** = 25th percentile of PD subject tremor power
 - **MID_THRESH** = 60th percentile
 - **HIGH_THRESH** = 85th percentile
 - Set **REG_THRESH** = median regularity score across PD windows
 - These 4 values are your classifier's only parameters — derived from data, explainable, adjustable
-

Stage 4 — Per-Window Classification

Rule-based classifier using the 4 thresholds:

Condition	Label
Power < LOW or dominant freq outside 4–6 Hz	None
Power ≥ LOW but < MID, or regularity < REG	Mild
Power ≥ MID but < HIGH	Moderate
Power ≥ HIGH and regularity ≥ REG	Severe

Each window gets one label. No probabilities needed for MVP.

Stage 5 — Hourly Aggregation

- Group windows by hour using session timestamps
- **Majority vote** across all windows in that hour → one bin per hour
- Track window count per hour:
 - ≥ 6 windows → confident
 - 3–5 windows → low confidence flag ⚠
 - < 3 windows → mark as no data, don't show a label

- Hours with no sessions at all → grey / no data in the view

Stage 6 — ON/OFF Derivation

Simple mapping on top of tremor bins — no separate classifier needed:

Tremor Bin	Motor State
------------	-------------

None	ON
------	----

Mild	ON
------	----

Moderate	OFF
----------	-----

Severe	OFF
--------	-----

This gives you the ON/OFF timeline as a byproduct of the tremor classifier, not a separate model.

Stage 7 — Validation Against UPDRS

Before building the view, sanity check the classifier:

- For SData subjects with UPDRS tremor subscores, check that:
 - High UPDRS tremor score → majority of their windows labeled Moderate/Severe
 - Low UPDRS score → majority labeled None/Mild
- This is not a formal accuracy metric — it's a face validity check for the hackathon
- Report the agreement rate as your "clinical grounding" claim in the pitch

Stage 8 — Output Structure

What gets passed to the clinician report view:

Per hour:

- hour timestamp

- tremor bin: None / Mild / Moderate / Severe
- motor state: ON / OFF
- window count + confidence flag
- dominant frequency (for tooltip detail)

Per day summary:

- % hours in each bin
- % ON-time
- peak tremor hour
- data coverage (hours with data / 24)

Stage 9 — Clinician Report View

7-Day Timeline Row (one row per day, 24 columns)

Color	Meaning
-------	---------

Green	ON (None/Mild tremor)
-------	-----------------------

Amber	Moderate
-------	----------

Red	Severe / OFF
-----	--------------

Grey	No data
------	---------

Daily Summary Card

- ON-time %
- Tremor severity breakdown as a small bar (% hours in each bin)
- Peak tremor period (e.g. "Most severe: 8–10am")
- Data coverage caveat if < 50% of hours have data

LLM Narrative Input

- Pass the daily summaries as structured JSON
 - The LLM converts them into the 3-paragraph clinical narrative
 - Example input signal: "Monday: 60% ON, peak Moderate tremor 8–10am, 6 hours of data"
-

What This Is Not Claiming

Important for pitch honesty:

- This is not a diagnostic tool — PD diagnosis is pre-established
- Thresholds are population-derived from SData, not clinically validated cutoffs
- Phase 2 replaces the rule-based thresholds with a trained classifier using UPDRS as ground truth
- Hourly bins assume session data is representative of that hour — not guaranteed with SData's call-based recordings